



**STERILE**  
**STERILE**  
**ALCOHOL**



### QUALITY STANDARDS

FG Clean Wipes is a pioneer in the manufacturing of wipes for cleanrooms. Quality is a prerequisite to everything we make. That means the dependability and performance of our sterile alcohol spray are unmatched.

### APPLICATIONS

- Surface cleaning
- USP <797> aseptic areas
- Use on gloved hands in sterile areas
- Critical medical device, biotechnology and pharmaceutical manufacturing areas
- Wipe downs for pass through controlled areas
- Ideally suited for use with FG Clean Wipes sterile dry wipes
- Cleaning of surfaces, such as counter tops, pass-throughs, and hoods

### BENEFITS

- Eliminates mixing and processing in-house
- Quick evaporation, leaving extremely low chemical residue
- Fully assembled, ready to use
- Economical pricing with bulk jugs
- Spray bottles provide custom engineered performance and safety
- Lot traceable

### FEATURES

- 70% USP Grade Isopropyl Alcohol
- 30% USP Purified DI Water or Water for Injection (WFI) Quality
- Alcohol blend filtered to 0.22 microns
- Each bottle is double bagged
- Pyrogen Free
- Gamma irradiated and validated to a sterility level of 10<sup>-6</sup>
- Certificate of Sterility provided with each case

### STERILE ALCOHOL ATTRIBUTES

<b>PH</b> 6.0 to 8.0	<b>APPEARANCE</b> Clear, colorless liquid
<b>ODOR</b> Characteristic IPA	<b>WATER QUALITY</b> USP Purified DI Water or WFI Quality
<b>DENSITY</b> 0.872 to 0.883 g/ml at 20°C	<b>PACKAGING</b> 16 fl. oz. spray bottle   12 bottles/case 32 fl. oz. spray bottle   12 bottles/case 1 gal. (3.7585L) bottle   4 bottles/case

ITEM #	MATERIALS	BOTTLE SIZE	PACKAGING
6-LS7030VSIPA-16B	70% IPA / 30% DI Water	16 oz Spray Bottle	12 bottles/case
6-LS7030VSWFI-16B	70% IPA / 30% WFIQ Water	16 oz Spray Bottle	12 bottles/case
6-LS7030VSIPA-32B	70% IPA / 30% DI Water	32 oz Spray Bottle	12 bottles/case
6-LS7030VSWFI-32B	70% IPA / 30% WFIQ Water	32 oz Spray Bottle	12 bottles/case
6-LS7030VSIPA-1G	70% IPA / 30% DI Water	1 gal Bottle	4 bottles/case



No matter your industry, we have the optimal critical cleaning solution.  
Find out more information at [www.FGCleanWipes.com](http://www.FGCleanWipes.com)

Phone: +1 (413) 598-8344  
Fax: +1 (413) 592-4112